

Pre-installation Checklist

Conditions for installation:

1. The ambient temperature:5°C - 35°C, and the relative humidity: $\leq 80\%$.
2. The laboratory should be free of violent vibration, airflow, strong electromagnetism and corrosive gases.
3. Equip the laboratory with proper work surface: 700mm (W) \times 2500 mm(L) \times 600mm(H).
4. Equip the laboratory with stable power supply of 220 V \pm 22 V / 50 \pm 1 Hz (grounded well). Two groups of power supply in different phase: Controlling power \geq 1 kW, and heating power \leq 4 kW. Reserve locations for one air switch and one three-core socket (20A,220V).

Preparation:

1. Quartz sand or fine sand (incinerated at 900 °C for half an hour), for measuring the petroleum products.
2. Absorbent cotton (medical), high temperature cotton (glass wool).

3. Oxygen: Prepared by industrial method (purity>99.5%).Do not use electrolytic oxygen. The cylinder pressure of gas supply should be more than 3MPa.

4. Please confirm connectors of the cylinder are consistent with Chinese standard (G5/8"-RHF) (the screw thread is on the outside) to make it match with the reducing valve the instrument equipped with, if not, please prepare the reducing valve (Guage for cylinder is 0-25MPa, guage for outlet is 0-1MPa) by yourself.

5. An electric balance with a resolution of 0.1 mg.

6. One desiccator.